



Experimental campaign on VEGA-3 laser: thin film targets investigation &

magnetic spectrometer testing F. Gatti^a*, F. Mirani^a, A. Maffini^a, M.S. Galli De Magistris^a, M. Galbiati^a, K. Ambrogioni^a,

D. Vavassori^a, D. Orecchia^a, D. Rastelli^b, A. Pola^a, J.L. Henares^c, A. Morabito^c, and M. Passoni^a

^aDepartment of Energy, Politecnico di Milano, Milano, Italy

- ^bRaylab s.r.l., Caravaggio (Bergamo), Italy
- ^cCentro de Láseres Pulsados, Salamanca, Spain

*francescol.gatti@polimi.it

Motivations and goals

Ion acceleration

Short laser pulses (10s fs) with relativistic intensities $(10^{18}-10^{22} \text{ W/cm}^2)$


$$a_0 = \frac{eE_0}{m_e\omega_0c} > 1$$

- Interaction with micrometric solid foils: Target Norma Sheat Acceleration (TNSA)
- **Plasma** generation: overdense regime $n_e > n_c$

$$n_c = \frac{m_e \omega_0^2 \epsilon_0}{e^2}$$

- Acceleration of high energy e⁻ and ions: 1-10s MeV [1]
- Enhanced TNSA [2]: addition of near-critical material $n_e \approx n_c$
 - Double Layer Targets (DLTs): few μ m solid + ~10 μ m low density -
 - > Volumetric laser-plasma interaction
 - Density: **few mg/cm**³

- Control and improve energy and number of particles [3]
- Potential compact and flexible sources
 - Application in **materials science** (e.g. PIXE [4, 5]), radioisotopes production, ...



dN/dE/dΩ


400

≥ 300 -

eg 200 -

100

 10^{10}




8 10 12 14 16 18 20

Proton energy [MeV]

Experimental setup

- VEGA-3 laser:
 - 30-50 fs
 - 12 μm FWHM
 - 22 J (4 J in spot on target)
 - ~1.25x10²⁰ W/cm²



Targets characterization soap layer ■ **Thickness** uniformity of Al targets assessed through EDXS maps in the range 750 nm - 6 µm Commercial Al Deposited Al **5%** max deviation 25% max deviation & up to 40% difference form nom. Avg. density $\longrightarrow > \sim 6.4 \text{ mg/cm}^3$ **Thickness** —> 9, 18, 36, 54 μ m C-foam Described in [6] Magnetron Sputtering (HiPIMS) deposition Deposited 1.5 μm Commercial 1.5 μm Deposited 0.75 μm Commercial 0.75 μm 3 μm


Carbon foam [2] on **1.5 μm and 6 μm Al**: Commercial 6 µm

Spectrometer design 1) Pinhole to select known small angle (3) Filter to remove heavy ions (2) Magnet to spreads ions according to energy 4) Detector arrays to collect signals Ion filter ■ Total size < 50 cm, with 5 cm, 0.4 T Energy range and resolution: **0.40±0.01** - 72±15 MeV Output signals (Ref. value: **10.2±0.8** MeV) to oscilloscope ■ Detectors' **pixels** ~2 mm² Charge Collection Efficiency modelled Voltage pulses amplitude readout — ▶ Real time Pixelated detectors (4) ■ Response matrix from GEANT4 **──** Absolute calibration Ion bunch Dipole magnet Inlet and pinhole 0.0 25 50 75 100 125 150 175 2

Target norma

Ion acceleration results

Deposited vs commercial



- Same E_{max} trend with thickness, peaked at ~3 µm
- Very similar cutoff energy
- Higher shot-to-shot reproducibility
- About **30% more** accelerated protons per unit solid angle

Single vs double layer Commercial, 6 µm solid Al ■ Trend of E_{max} 54 µm foam independent of substrate ■ Peak between 18-36 18 μm foam μm foam thickness About 2 MeV higher 0 μm foam cutoff energy 10^{9} ■ About **70% more** accelerated protons per unit solid angle E [MeV] ▲ Commercial, TP (6 μm) Δ Deposited, TP (1.5 μm) /Su 150 - G O 100 - 2

20 40 60

Foam Th. [µm]

Measurements and comparisons 10^{13} ■ **First test** in laser-driven acceleration experiment **└** Spectrometer Front-end electronics sized independently for RCF each channel [1/MeV/sr]

- Very good qualitative agreement with Thomson
- Parabola spectra
- Quantitative agreement with RCFs data
- Cross-calibration of Thomson Parabola ■ Test without magnetic field: no unwanted signals recorded (e.g. EMPs, _{0.5 MeV} photons) C-ions B-field \odot

erc









Conclusions and perspective

20

Foam Th. [μm]

40

- Commercial rolled targets suffer of strong fluctuation of thickness and surface quality, with detrimental effects on accelerated ions
- Optimized HiPIMS-deposited targets can provide better shot-to**shot stability** of accelerated protons
- Carbon foams of optimal thickness enhanced both max energy and particle number, while the solild layer is still of much relevance
- The developed spectrometer provided spectra in accordance with the other diagnostics employed
- The retrieved proton numbers were in accordance with literature results and were also used to cross-calibrate the Thomson Parabola
- Improvement in spectrometer

E [MeV]

 $B = 0.4 \text{ T} -- 1^{st} \text{ exp}$ $B = 0.0 \text{ T} \dots 2^{nd} \text{ exp}$

0.0 0.2 0.4 0.6 0.8 1.0

- electronics and data
- **Experimental campaign** @ **e**i scheduled for October

 \odot

- Materials characterization with protons (PIXE)
- Use of both deposited targets and spectrometer

Bibliography:

Modulated stopper deposited with HiPIMS

[1] A. Macchi, et al., Rev. Mod. Phys. 85.2 (2013): 751.

Single photodiode

(~ 1 mm size)

- [2] A. Maffini, et al., Appl. Surf. Sci. 599 (2022): 153859. [3] I. Prencipe, et al., NeW J. Phys. 23.9 (2021): 093015.
- [4] M. Passoni, et al., Sci. Rep. 9.1 (2019): 9202.
- [5] F. Mirani, et al., Sci. Adv. 7.3 (2021): eabc8660.
- [6] A. Maffini, et al., EPJ tech. instrum. 10.1 (2023): 15.
- [7] F. Gatti, et al., IEEE Trans. Instrum. Meas., Accepted.